

Antibody Staining Protocol with TSA amplification

- 1) Wash embryos 2x 30 minutes in PBTw with rocking.
- 2) Block for 30 minutes in 100 μ L PBTw + 5% NGS with rocking.
- 3) Add the primary antibody (at the recommended dilution in PBTw +NGS). Make sure the embryos are properly mixed.
- 4) Incubate overnight at 4°C.
- 5) Recover the primary antibody to use again by transferring to another tube.
- 6) Wash embryos 3x quickly (let embryos sink, then remove liquid) with PBTw.
- 7) Wash 3x 30 minutes with PBTw with rocking.
- 8) Block 30 minutes in 100 μ L PBTw + 5% NGS with rocking.
- 9) Add biotinylated secondary antibody (at the recommended dilution in PBTw +NGS, usually 1:300).
- 10) Incubate at room temperature for 60 minutes (no rocking necessary).
- 11) Wash 3x quick with PBTw.
- 12) Wash 3x 30 minutes with PBTw.
- 13) Block for 30 minutes in 100 μ L PBTw + 5% NGS with rocking.
- 14) Incubate 60 minutes in streptavidin-HRP (Jackson Labs 1:100 dilution).
- 15) Wash 3x quick with PTw.
- 16) Wash 3x 30 minutes with PBTw. BE SURE TO REMOVE ALL PTw AFTER LAST WASH.
- 17) Incubate 15-20 minutes in TSA reagent (TSA-Cy3 or TSA-Cy5 diluted 1:50 in TSA amplification diluent).
- 18) Wash 3x quick with PBTw.
- 19) Wash 3x 10 minutes with PBTw.
- 20) Embryos are ready for the next step (ie mounting or the next antibody).

Solutions:

PBTw:

- 1X PBS
- 0.1% Tween-20
- 0.1% BSA (bovine serum albumin)

PBTw + NGS:

- 1X PBS
- 0.1% Tween-20
- 0.1% BSA (bovine serum albumin)
- 5.0% NGS (normal goat serum)

TSA kits (Perkin Elmer)

- NEL-705A – TSA Cy5 kit
- NEL-704A – TSA Cy3 kit
- NEL-701 – TSA FITC kit