

In vitro transcription protocol

1. PCR with appropriate primer set to obtain cDNA fragment flanked by RNA polymerase binding sites.

We generally run 25 ul PCR reactions.

2. Check PCR product size and quantity on agarose gel.

Run 2.5 ul on a normal agarose gel

3. Precipitate PCR product.

20 ul of PCR product (see note)
 10 ul 3M Sodium Acetate (pH5.5)
 70 ul water
 900 ul cold 100% ethanol

- incubate at -20C for 30 minutes and spin at max speed for 30 minutes
- remove ethanol from RNA pellet and wash pellet with 70% ethanol (500ul).
- resuspend in 10 ul of DEPC water

Note: the goal is to have 1ug of DNA template, however in practice this much is rarely needed to make a good probe. 20 ul of a DNA sample that shows a decent intensity band on a gel will be fine. You may have to use your own judgment.

4a. In vitro transcription reaction T7 or SP6 (we use NEB enzymes):

PCR product	10 ul	(use 1 ug of DNA in DEPC dH2O)
10X transcription buffer	2 ul	
10X DIG labeling mix	2 ul	(Roche, contains DIG-UTP)
RNasin (40U/ul)	0.5 ul	(Promega, prevents RNase activity)
RNA polymerase (15-17U/ul)	1 ul	(T7 or SP6)
DEPC dH ₂ O	<u>4.5 ul</u>	
	20 ul	

- incubate 2-3 hours at 37°C.
- add 1 ul of RNasin (40U/ul) and 1ul of RQ1 DNase.
- incubate 30 minutes at 37C. can store at -80°C.

4b. In vitro transcription reaction T3 (we use Promega enzyme):

PCR product	10 ul	(use 1 ug of DNA in DEPC dH2O)
5X transcription buffer	4 ul	(Promega)
0.1M DTT	2 ul	
10X DIG labeling mix	2 ul	(Roche, contains DIG-UTP)
RNasin (40U/ul)	0.5 ul	(prevents RNase activity)
RNA polymerase (15-17U/ul)	1 ul	(T3)
DEPC dH ₂ O	<u>0.5 ul</u>	
	20 ul	

- incubate 2-3 hours at 37°C.
- add 1 ul of RNasin (40U/ul) and 1ul of RQ1 DNase.
- incubate 30 minutes at 37°C. can store at -80°C.

Note, the protocol is the same for biotinylated probes. Just substitute Roche biotin RNA labeling mix in the above protocol.

5. Check RNA agarose gel (1%)

We normally run 1-2 ul of the RNA reaction on a gel to examine 1) if RNA was made and 2) to get a rough idea of the concentration. Running a formaldehyde gel is optional, regular gels work just fine for this purpose. Normally, we run the gel while performing step 6. However, the RNA reaction can be stored at -80C while you run the gel and then you can proceed to step 6 once you are confident about the probe reaction working.

6. Probe fragmentation and precipitation.

Some protocols suggest that the optimal fragment size for RNA hybridization is 500bp. We have found that hydrolyzed probes work slightly better than non-hydrolyzed probes. If your probe is <1kb bp is not a good idea to hydrolyze as your fragments will be very small. Hydrolysis of larger probes (>1kb) can be achieved by using a mild alkaline solution as follows:

- add an equal volume (20 ul) of 0.2M Na₂CO₃ pH 10.2 to each probe.
- incubate at 60°C for 15 minutes
- Place on ice, add 20ul of 7.5M NH₄OAc and 160 ul of 100% ethanol.
- incubate at RT for 10 minutes and spin at max speed for 30 minutes
- remove ethanol from RNA pellet and wash pellet with 70% ethanol (500ul).
- resuspend in 150 ul of resuspension buffer and store at -80°C.

Solutions:

3M Sodium Acetate (pH 5.5)

100% and 70% Ethanol

0.2M Na₂CO₃ (pH 10.2)

7.5M Ammonium Acetate

Resuspension buffer: 50% Formamide, 5mM Tris (pH7.5), 0.5mM EDTA, 0.01% Tween

Roche Dig RNA labeling Mix

T7, T3, and SP6 enzymes and accompanying buffers